# TekniZorb™ LXS Meltblown Polypropylene

TekniZorb<sup>™</sup> LXS offers a very low-shedding wiping solution to deposit IPA while effectively removing residue from the surface in your controlled environment.

# APPLICATIONS

- Ideal for efficient surface cleaning and sanitization between and during processing.
- Useful for pass through wiping and PPE wipe down in gowning area.
- Suitable for use in ISO 5 and higher cleanrooms

# ■ FEATURES & BENEFITS

- Excellent chemical compatibility
- Applies IPA efficiently while minimizing VOC emissions
- Exceptionally clean
- Available in resealable pouches



Teknipure

# TekniZorb™ LXS Meltblown Polypropylene Wipers

Sterile and Non-Sterile

		Edge				Wiper	
Product	Class	Finish	<b>Pre-Sat</b>	Sterile <sup>†</sup>	IPA%*	Size	Packaging
TS1MPI70-911	ISO 5	Standard	<b>•</b>		70%	9" x 11" (23cm x 28cm)	50 ea/pouch; 24 pouches/case
TS1MPI70-911S	ISO 5	Standard	<b>•</b>	<b>♦</b>	70%	9" x 11" (23cm x 28cm)	50 sheets/pouch; 24 pouches/case

<sup>†</sup>Validated sterilization via gamma irradiation to a Sterility Assurance Level of 10<sup>-6</sup> according to AAMI Guidelines

## **Physical Properties**

<b>Basis Weight</b>	35 g/m <sup>2</sup>	
Material	100% Polypropylene	

### **Purity Specifications**

Particles & Fibers	Target (x 10 <sup>6</sup> particles/m <sup>2</sup> )
Particles >0.50um	<20.0
Particles >5.0um	<3.0
Fibers >100um	<1000

<sup>\*</sup>Particles: IEST-RP-CC004.3 Section 6.1.4

Fibers: IEST-RP-CC004.3 Section 6.2.2 (Orbital Shake Test)

#### **Shelf Life**

Non-sterile Dry Wipes: 5-years from the date of manufacture. Non-sterile Pre-wetted Wipes: 3-years from the date of manufacture. Sterile Dry & Pre-wetted Wipes: 3-years from the date of manufacture.

#### Nonvolatile Residue

Extractant	Target
IPA	<0.005 g/m <sup>2</sup>
DI Water	<0.010 g/m <sup>2</sup>

<sup>\*</sup>IEST-RP-CC004.3 Section 7.1.2 (Short Term Extraction)

#### **Extractable Ions**

lon	Target		
Sodium (Na+)	<0.50 ppm		
Potassium(K+)	<0.50 ppm		
Calcium (Ca <sup>2+</sup> )	<1.0 ppm		
Chloride (Cl <sup>-</sup> )	<1.0 ppm		
Magnesium (Mg <sup>2+</sup> )	<1.0 ppm		
* 1505 55 6000 1 5 600 15 (6)			

<sup>\*</sup> IEST-RP-CC004.3 Section 7.2.2.1B (Standard Extraction)





Toll Free: 844.309.2376





<sup>\*</sup>IPA (Isopropyl Alcohol) meets or exceeds USP purity requirements